PHYSICIAN MANUAL	ALPHA (S) & ALPHA (M) BONE CONDUCTION HEARING SYSTEMS
MANUAL DEL MÉDICO	ALPHA (S) Y ALPHA (M) SISTEMAS AUDITIVOS POR CONDUCCIÓN ÓSEA
MANUEL DU MÉDECIN	ALPHA (S) ET ALPHA (M) SYSTÈME AUDITIF À CONDUCTION OSSEUSE
HANDBUCH FÜR DEN ARZT	ALPHA (S) & ALPHA (M) KNOCHENLEITUNGSHÖRSYSTEME
MANUALE PER IL MEDICO	ALPHA (S) E ALPHA (M) PROTESI ACUSTICHE A CONDUZIONE OSSEA
HANDLEIDING VOOR ARTSEN	ALPHA (S) & ALPHA (M) BEENGELEIDINGSHOORSYSTEMEN

CE 0086

Federal Law (USA) restricts this device to sale by or on the order of a practitioner licensed by the law of the State in which he/she practices to use or order the use of the device.

1 ABOUT THIS MANUAL

This Physician Manual provides a description of the Alpha (S) & Alpha (M) Bone Conduction Hearing Systems. This manual includes a description of the system's features, performance, intended use, and important safety information. This manual provides the information necessary for the physician to prescribe the Alpha (S) & Alpha (M) Bone Conduction Hearing Systems to the appropriate patient population. Please read this manual thoroughly before prescribing the system.

A Patient Instruction Manual is included with each Alpha Sound Processor and should be referred to for information regarding:

- Kit Contents
- Instructions for Use
- Patient Warnings and Precautions
- Important Things to Remember
- Maintenance
- Customer Support

2 INTRODUCTION

Warnings and Cautions Overview

The following section contains a description of all physician related Warnings and Cautions associated with the Otomag™ Bone Conduction Hearing Systems. Additional patient and use related Warnings and Cautions can be found in the Patient User Manual. Knowledge of the following information is essential for safe operation of the device. Specific Warnings and Cautions pertaining to physician instructions are included in the body of the manual.

Warnings

- Thoroughly read the Physician Manual and the Patient Manual before prescribing the Otomag™ System.
- To avoid potential performance and / or safety issues, adhere carefully to patient selection criteria
- A blow to the head in the area of the Otomag™ Alpha Sound Processor may result in damage to the device or device failure.
- When using the Otomag™ Alpha (S) Softband, make sure that a finger can be placed between the skull and Softband to ensure that the Softband is not too tight.
- Do not perform Magnetic Resonance Imaging (MRI) on a patient wearing an Otomag™ Alpha Sound Processor, Magnetic Spacer, Headband or Softband. The external portions of the device must be removed prior to entering the MR environment. Severe injury may result if the external portions of the device are not removed prior to entering the MR environment.
- The Alpha (M) Magnetic Implant can be scanned safely using MRI only under certain conditions.
- Carefully consider the potential patient's medical history to rule out allergies and other reactions to device materials prior to prescribing the Otomag™ System.
- The Alpha (M) Magnetic Implant, the Implant Template and Surgical Template are supplied on a plastic card provided sterile in gas-permeable packaging. Handle the package with care. Severe impact may rupture the sterile package.
- The implant sterile package should be stored at normal room temperature. The package may be stored at temperatures between -4 °F and +122 °F (-20 °C and +50 °C).
- The 'use by' date is stamped on the sterile package. If it has expired, return the device to Sophono.
- The Alpha (M) Magnetic Implant, the Implant Template and Surgical Template are single use items. The sterile package contains information indicating ethylene oxide processing. Before opening the sterile package, inspect it carefully. If the package is ruptured, or exposure to ethylene oxide processing is not indicated, return the package to Sophono.
- The Magnetic Implant and Implant Template are re-sterilizable according to the provided instructions. Failure to sterilize appropriately can result in serious infection.
- The Surgical Template is *not* re-sterilizable.
- If re-sterilization is necessary, the Alpha (M) Magnetic Implant and Implant Template are not intended to be flash steam sterilized. Flash steam sterilization has not been validated to effectively sterilize the components.
- Do not re-sterilize the Alpha (M) Magnetic Implant more than a total of three (3) times. Additional sterilization cycles may damage the implant.
- Purchased screws may be provided non-sterile and if so must be sterilized prior to implantation according to the manufacturer's instructions. Failure to sterilize appropriately can result in serious infection.
- Excessive, long term use of magnetic spacers that are too strong may lead to pain and or pressure necrosis for the patient. Always recommend to the patient the lowest magnetic strength spacer that provides sufficient retention.
- Do not use the Alpha (M) configuration of the device if the skin contact area of the sound processor is compromised in any way (i.e. – unhealed incision, psoriasis, or eczema). Failure to follow this warning may result in prolonged healing, further irritation, or infection.
- Both the Alpha (S) and Alpha (M) configurations of the device contain magnets which could impact the performance of magnetically sensitive medical devices such as a pacemaker, an implantable cardioverter defibrillator (ICD), or ventricular shunt. The magnets of the softband/headband of the Alpha (S) and the magnetic spacer of the Alpha (M) should therefore be kept away from such devices.

Precautions

- The Alpha Sound Processor, Magnetic Implant, Headband, and Softband are intended for single-patient use only. Do not use any of these products if they have previously been placed in or with another patient. Reuse of the magnetic implant carries a serious risk of severe infection and other major medical complications.
- The patient should be instructed to select a lower strength spacer if the pressure on the skin is uncomfortable or pain, numbness, or tingling occurs. If discomfort continues, advise the patient to discontinue wear and visit the clinic for evaluation and refitting.
- Patients with a pacemaker or ICD should consult their care provider for safe use of the Otomag system with their cardiac medical device.
- Keep Otomag™ System components away from magnetic data storage media and electronic devices. The Alpha Sound Processor, Magnetic Implant, Headband, and Softband all contain magnets that can potentially damage magnetic data storage media and electronic devices.
- To avoid damage to the device, do not subject the Alpha Sound Processor to temperatures below -4°F (-20°C) or above 122°F (+50°C) at any time.
- To avoid potential damage to the device, do not store the Otomag™ System near magnets or magnetic devices other than the magnetic spacer and implant.
- Prior to use, inspect packaging to verify that it is intact and has not been damaged. If the packaging is damaged, do not use the device and return it to Sophono. Please see the address and phone number for Sophono at the end of this Physician Manual.
- Refer to the original manufacturer's instructions for use for specific directions on how to use the implantable screws.
- Refer to the original manufacturer's instructions for use for specific directions on how to appropriately use the sterilization pouches.

Product Description

The Otomag™ Bone Conduction Hearing System is a family of sound processors and accessories that operate on the principle of bone conduction of sound vibrations.

The Otomag™ System is configured in either of two configurations. The first configuration is the Alpha (S), where the Otomag™ Sound Processor is attached magnetically to a Headband or Softband. The second configuration is Alpha (M), where the Otomag™ Sound Processor is attached magnetically to an implanted magnet. The Headband, Softband, or Magnetic Implant holds the sound processor against the head, and vibration is transduced through direct contact with the patient's skin and the bone below.

The Otomag™ System is designed for use for those patients with conductive hearing loss, those patients who have sensorineural hearing loss up to 45 dB in combination with their conductive loss, and single sided deafness as defined in the indications for use. The prescriptive formula and adjustments available to the audiologist in the software allow for programming the Otomag™ System for individual patient hearing loss.

Indications for Use

The Otomag™ Alpha Sound Processor is intended for use with the Otomag™ Headband or Otomag™ Softband (no age limitations), or with the Otomag™ Magnetic Implant (patients 5 years of age and up) for the following patients and indications:

- Patients with conductive or mixed hearing losses, who can still benefit from amplification of sound. The pure tone average (PTA) bone conduction (BC) threshold for the indicated ear should be better than 45 dB HL (measured at 0.5, 1, 2, and 3 kHz).
- Bilateral fitting is applicable for most patients having a symmetrically conductive or mixed hearing loss. The difference between the left and right sides' BC thresholds should be less than 10dB on average measured at 0.5, 1, 2, and 4 kHz, or less than 15 dB at individual frequencies.
- Patients who have a profound sensorineural hearing loss in one ear and normal hearing in the opposite ear, who for some reason will not or cannot use an AC CROS. The pure tone average (PTA) air conduction (AC) threshold of the hearing ear should be better than 20 dB HL (measured at 0.5, 1, 2 and 3 kHz).

Contraindications

Contraindications for the Otomag™ Alpha Sound Processor connected to the Otomag™ Headband, Softband, or Magnetic Implant are as follows:

Any factor that would cause a clinician to refer the patient for medical assessment will temporarily, or in some cases permanently, halt the process of hearing aid fitting. These factors include:

- A hearing loss of sudden onset;
- A rapidly progressing hearing loss;
- Pain in either ear;
- Tinnitus of sudden recent onset, or unilateral tinnitus;
- Unilateral or markedly asymmetrical hearing loss of unknown origin;
- Vertigo (e.g. dizziness)

Patient Selection Criteria

- Patients to be fitted with the Otomag™ Alpha (S) Sound Processor and its Headband or Softband may be any age.
- Patients to be fitted with the Otomag™ Alpha (M) Sound Processor and its Magnetic Implant should be 5 years old or older.
- Patients and/or patient guardians must be able to understand that the Otomag™ System will not restore hearing to normal and will not prevent additional hearing impairment resulting from conditions such as disease, old age hearing loss, genetic hearing loss, ototoxic drugs or noise induced hearing loss.

Possible Adverse Effects**General**

It is possible that a reaction to the materials of the Headband, Softband, or Magnetic Implant could occur. Such a circumstance may necessitate that use of the device be discontinued. Please report any adverse effects to Sophono. The address and phone number for Sophono are listed at the end of this Physician Manual. If a reaction to the Headband or Softband occurs, Sophono will supply appropriate instructions so that the device can be sent to the laboratory for analysis.

Risks that may occur during surgery:

The risks involved in surgery and complications related to surgery include the normal risks of the operative procedure and general anesthesia. These include pain, blood clots and/or circulatory problems, infection, cardiac arrest, and even death.

Warning

- To avoid potential performance and / or safety issues, adhere carefully to patient selection criteria.
- Carefully consider the potential patient's medical history to rule out allergies and other reactions to device materials prior to prescribing the Otomag™ System.

Risks that may occur after surgery:

Potential risks after surgery include inflammation and complications associated with general anesthesia. Potential risks related to the implant following surgery include infection, reparative granuloma and inflammatory reaction.

Specific risks related to operating on the mastoid include infection, inflammation, numbness or stiffness about the ear, dural exposure, Cerebral Spinal Fluid (CSF) leak, perforation of sigmoid sinus, subdural hematoma, disturbance of taste or balance, and noticeable changes in head noise.

It is possible that a reaction to the implant materials and or failure of the implanted device can occur and may necessitate an operation to remove and/or replace the implant.

Warning

- A blow to the head in the area of the Otomag™ Alpha Sound Processor may result in damage to the device or device failure.
- When using the Otomag™ Alpha(S) Softband, make sure that a finger can be placed between the skull and Softband to ensure that the Softband is not too tight.
- Do not perform MRI on a patient wearing an Otomag™ Alpha Sound Processor, Magnetic Spacer, Headband or Softband. The external portions of the device must be removed prior to entering the MR environment. Severe injury may result if the external portions of the device are not removed prior to entering the MR environment.
- The Alpha (M) Magnetic Implant can be scanned safely using MRI only under certain conditions.

MRI Information

Non-clinical testing demonstrated the Alpha (M) Magnetic Implant is MR Conditional and can be scanned safely using MRI only under the following conditions:

- remove all external components including the Otomag Alpha Sound Processor, Magnetic Spacer, Headband or Softband before entering the MR environment
- static magnetic field of 3 Tesla or less
- spatial gradient field of 720 Gauss/cm or less
- maximum whole-body averaged specific absorption rate (SAR) of 4 W/kg in the First Level Controlled Mode for a maximum scan time of 15 minutes of continuous scanning.

In non-clinical testing the Alpha (M) Magnetic Implant produced a maximum temperature rise less than 3.2°C during 15 minutes of continuous MR scanning in the First Level Controlled Mode at a maximum whole-body averaged SAR of 4 W/kg.

The computed implant temperature increase in response to the worst-case, time-averaged gradient field (94.7 Tesla/sec) for the 15 minute possible exposure for a series of clinical MRI scans is less than 2.6°C.

Image Artifact

The maximum image artifact size extends approximately 5-cm relative to the size and shape of the implant when scanned in nonclinical testing using the Gradient echo (GRE) pulse sequence in a 3 Tesla/128-MH, Excite, Software 14X.M5, General Electric Healthcare, Milwaukee, WI; active-shielded, horizontal field MR system with a send receive RFcoil. Artifacts were less with the T1-weighted, spin echo pulse sequence.

Implant Function Following MR Scanning

In non-clinical testing, when positioned parallel to the patient table, the Alpha (M) Magnetic Implant maintained over 95% of its original magnetic strength after 10 insertions into a static MRI field and over 10 minutes of pulse sequence in a 3 Tesla Siemens Tri Clinical MRI Scanner (MRC20587). To minimize demagnetization of the internal magnets, Sophono recommends that the patient's head is positioned along the long axis of the MRI scanner and that the patient is instructed to avoid head movement while lying on the patient table. If the head is tilted at an angle to the long axis of the patient table during a scan, the internal magnets may be demagnetized and a stronger external magnet may be needed to keep the device in place.

Caution

- After MRI exposure the Magnetic Spacer Strength should be adjusted as necessary following the instructions "Instructions for Placement / Sizing of the Magnetic Spacers and Alpha Sound Processor" included in this manual.

Symbols on Equipment and Labeling

Symbol	Description
	Consult Physician Manual
	Batch Code or Lot Number
	Serial Number
	Keep Away From Rain / Keep Dry
	Manufacturer
	MR Unsafe: The Alpha Sound Processor, Magnetic Spacer, Headband or Softband must be removed before entering the MR environment
	MR Conditional: The Alpha (M) Magnetic Implant can be scanned safely using MRI only under the following conditions: <ul style="list-style-type: none"> remove all external components including the Otomag Alpha Sound Processor, Magnetic Spacer, Headband or Softband before entering the MR environment static magnetic field of 3 Tesla or less spatial gradient field of 720 Gauss/cm or less maximum whole-body averaged specific absorption rate (SAR) of 4 W/kg in the First Level Controlled Mode for a maximum scan time of 15 minutes of continuous scanning.
	Product is Latex Free
	Interference may occur in the area of equipment which is labeled with the following symbol.
	Temperature Limitation: Storage and/or use limited to between -4°F (-20°C) and 122°F (+50°C)
	Sterile by Ethylene Oxide. The Alpha (M) Magnetic Implant, the Implant Template and Surgical Template are provided sterile.
	Single Use
	Authorized representative in the European Community

3 PATIENT FITTING PROCEDURES

Alpha (S) System: Instructions for Placement of the Softband or Headband and Alpha Sound Processor

Notes

- Record the product and lot or serial numbers for the Alpha Sound Processor in the patient file and send the patient to a qualified audiologist for proper programming of the Alpha Sound Processor.

The Otomag™ Softband is an elastic Softband similar to a sweat band used in sports. The Alpha (S) Sound Processor magnetically adheres to the Softband via a magnetic pad sewn into the band. No surgery is required and the Alpha (S) can be used immediately after being programmed by an audiologist.

- To use the Otomag™ Softband, attach the Alpha (S) Sound Processor to the magnetic pad on the Softband.
- Test that the Alpha works by putting it on yourself, plugging your ears, and introducing sound.
- Put the Softband around the patient's head, quite loosely at first by slipping the end of the band through the slot next to the magnetic pad and attaching the mating Velcro® loops and hooks to each other.
- Place the magnetic pad against the mastoid or other boney location.
- Tighten the Softband by adjusting the Velcro® strap until it is close-fitting enough to ensure effective sound transmission, but not tight enough to cause discomfort. The Softband is too loose if it slips off of the head, can be rotated easily around the head or no sound transmission can be heard. The Softband is too tight if the patient complains of pain or discomfort while wearing the Softband. Make sure that a finger can be placed between the skull and Softband to ensure that the Softband is not too tight.

Warning

- When using the Otomag™ Alpha (S) Softband, make sure that a finger can be placed between the skull and Softband to ensure that the Softband is not too tight.

The Otomag™ Headband is a springy Headband, similar to that used for earphone devices. It comes in one size and is not adjustable. The Alpha (S) Sound Processor magnetically adheres to the Headband via a magnetic pad fitted to the Headband. No surgery is required and the Alpha (S) can be used immediately after being programmed by the audiologist.

- To use the Otomag™ Headband, attach the Alpha (S) Sound Processor to the magnetic pad on the Headband.
- Test that the Alpha works by putting it on yourself, plugging your ears, and introducing sound.
- Put the Headband on the patient's head and place the magnetic pad against the mastoid or other boney location.

The choice of a Headband or Softband is a matter of patient preference.

Alpha (M) System: Compatible Screws

Notes

- The Alpha (M) implant requires 5 screws for implantation that are **NOT** included with the implant.

Screws

The Alpha (M) implant requires 5 screws to be used to affix the implant to the skull. These screws are NOT included as part of the Alpha (M) implant system and must be obtained separately.

The implant bracket holes are sized to accept Stryker/Leibinger 1.7 mm diameter Craniomaxillofacial titanium Self Tapping Screws, (Gold) and 1.7 mm diameter Self Drilling Screws (Silver) in 4 mm lengths and 1.9 mm Emergency Screws (Purple) in 3 mm length.

Craniomaxillofacial Titanium Screws	1.7 mm dia. x 4.0 mm Self Tapping	1.7 mm dia. x 4.0 mm Self Drilling	1.9 mm dia. x 3.0 mm Emergency
Stryker/Leibinger Part Number	50-17004	50-17904	50-19003
U.S. FDA 510(k) Clearance Number	K970912	K970912	K963739

While all screw sizes / types specified above have been validated for this application, the Stryker/Leibinger 1.7mm x 4 mm self-drilling screw, part number 50-17904, is considered optimal.

Other commercially available craniomaxillofacial titanium screws may be used as long as the screw shaft diameter is less than 2.0 mm and the screw head size is greater than 2.45 mm. Screws must be 4 mm in length (3mm for emergency screws).

The screws must be sterile prior to implantation. If the screws are not purchased sterile, refer to the screw manufacturer's instructions for sterilization.

Warning

- Purchased screws may be provided non-sterile and if so must be sterilized prior to implantation according to the manufacturer's instructions. Failure to sterilize appropriately can result in serious infection.

Caution

- Refer to the screw's original manufacturer's instructions for use for specific directions on how to use the screws.

If the initial attempt to place a screw does not work and the screw hole becomes "stripped" an emergency "rescue" screw can be threaded into the stripped hole. The largest size screw that can be used is 1.9 mm in diameter.

Alpha (M) System: Information Regarding Sterilization, Storage, Re-Sterilization and Placement of the Magnetic Implant

Sterilization and Storage:

Warning

- The Alpha (M) Magnetic Implant, the Implant Template and Surgical Template are supplied on a plastic card provided sterile in gas-permeable packaging. Handle the package with care. Severe impact may rupture the sterile package.
- The sterile package should be stored at normal room temperature. The package may be stored at temperatures between -4 °F and +122 °F (-20 °C and +50 °C).
- The 'use by' date is stamped on the sterile package. If it has expired, return the device to Sophono.
- The Alpha (M) Magnetic Implant, the Implant Template and Surgical Template are single use items. The sterile package contains information indicating ethylene oxide processing. Before opening the sterile package, inspect it carefully. If the package is ruptured, or exposure to ethylene oxide processing is not indicated, return the package to Sophono.

Re-Sterilization Instructions:

Warning

- The Magnetic Implant and Implant Template are resterilizable according to the provided instructions. Failure to sterilize appropriately can result in serious infection.
- The Surgical Template is *not* re-sterilizable.
- If re-sterilization is necessary, the Alpha (M) Magnetic Implant and Implant Template are not intended to be flash steam sterilized. Flash steam sterilization has not been validated to effectively sterilize the components.
- Do not re-sterilize the Alpha (M) Magnetic Implant more than a total of three (3) times. Additional sterilization cycles may damage the implant.

The Alpha (M) Implant Template is an accessory to assist the surgeon in determining the correct diameter and depth of the bone bed. To prepare the accessory for steam sterilization, use the following steps:

- Transfer the template into the autoclave pouch.
- Place a single steam sterilization indicator into the pouch with the template.
- Immediately seal the autoclave pouch.

To prepare the Alpha (M) Magnetic Implant for steam sterilization, use the following steps:

- Transfer the implant into the autoclave pouch.
- Place a single steam sterilization indicator into the pouch with the implant.
- Immediately seal the autoclave pouch.

Exact compliance is required with the sterilization equipment manufacturer's instructions for use, maintenance, and periodic monitoring. The following gravity steam sterilization (autoclaving) parameters have been validated for sterilization of the implant and implant template.

- Limitations on Reprocessing: Magnetic implants may be re-sterilized three additional times if not used during the implant procedure. Implants should be stored in the sterilization pouch. However, the implant must be transferred to a new sterilization pouch just prior to re-sterilization and implantation.
- Cycle Type: Gravity
- Sterilization Time: 15 minutes minimum
- Sterilization Cycle Temperature: 132 degrees C / 270 degrees F
- Drying Time: 30 minutes minimum
- Packaging: Cardinal Health Self Sealing Autoclave Pouch P/N 92510 (or equivalent)
- Storage in Sterile State: The Magnetic Implant has not been validated for storage after sterilization.

Upon completion of sterilization of each component, verify that both sterilization indicator strips (interior and on exterior of package) indicate the steam cycle was completed correctly. Follow manufacturer instructions for interpretation of these indicators. If both the indicators do not show completion of the sterilization cycle, repeat the cycle.

Sterilization Pouch Specifications

To ensure the validated parameters are replicated appropriately, the following (or equivalent) is to be used to individually pouch each implant prior to performing the sterilization cycle:

- Manufacturer: Cardinal Health, <http://www.cardinal.com>, 800.964.5227
- Catalog Number: 92510
- Size: 5.5" x 10" Paper Self Seal Pouch
- Materials: Arjowiggins 60g/m2 Paper / CN PET 12/CP40 Micron
- Total Product Area: 0.035 m2
- Paper Area: 0.035 m2
- Ink: Medical Grade

Caution

- Refer to the original manufacturer's instructions for use for specific directions on how to appropriately use these sterilization pouches.

Internal Steam Indicator Strip

As described above, each time an implant is steam sterilized, place a single indicator strip labeled for use with the specified gravity steam sterilization cycle into the inside of the specified pouch. Follow manufacturer's instructions for interpretation of results.

Placement of the Magnetic Implant:

Pre-Surgical Training Requirements

Before attempting to implant the Otomag™ Implant, surgeons are required to complete training provided by Sophono. During the training they will receive guidance regarding the procedure, a thorough understanding of the device and surgical procedure and other important information. Contact Sophono Customer Services at the contact number listed in Section 6 of this manual to arrange training.

Surgical Risks

The implantation of the Otomag™ Implant requires a surgical procedure involving the removal of tissue and bone. As with any surgical procedure involving a general anesthetic, certain risks exist. Those risks include the risk of infection from poor hygiene or non-sterile technique, irritation of tissues during the implantation procedure, unintended cuts and abrasions, and irritation of tissues due to unintended exposure to implant materials. The implant is made of materials considered to be hypo-allergenic but in the case of unintended damage to the device, a patient could come in contact with substances that could generate an adverse reaction.

Surgical Site Preparation:

The hair should be clipped and the scalp shaved at the incision site immediately before surgery. With the incision site skin free of gross contamination, apply antiseptic (such as iodophors, alcohol-containing products, or chlorhexidine gluconate) in the area of the proposed incision to an area large enough to allow for the extension of the incision or creation of a new incision or drain sites (if necessary).

Pre-Surgical Antibiotic Prophylaxis:

Antibiotic prophylaxis should be considered if the surgeon determines that this surgical procedure may be associated with a considerable risk of infection or if a post-operative infection would cause a serious hazard to the patient's recovery or well being.

Surgical Implantation Steps

Step 1 - Making The Incision

The proper position for the implant is 6 cm from the ear canal at a 45-degree angle posterior and superior to the ear canal. The incision should be made 7.5 to 8 cm from the external ear canal (1.5 to 2 cm from the location of the implant so the incision is not directly over the implant). The Surgical Template will aid the surgeon in ensuring proper placement by marking the skin with a surgical marking pen. After making the incision, expose the skull where the implant will be placed. The incision must be large enough to allow placement of the implant. See Figure 1.


Figure 1: Incision Location

Figure 2: Bone Beds


Figure 3: Implant Installation Diagram

Step 2 - Create a bone bed (optional)

The implant is approximately 2.6 mm thick and optimal tissue thickness over the implant is 4mm to 5 mm. The Surgical Template is provided with 6 mm and 4 mm gauges for determining the appropriate skin flap thickness. Bone beds are recommended for the implant if the tissue thickness over the implant is less than 3 mm. See Figure 2.

To create a bone bed, mark the position of the implant magnets on the bone using the Surgical Template and drill the bone bed to a depth of approximately 3 mm. The implant magnets have a diameter of 10 mm. The Implant Template is useful in ensuring that the bone beds are of sufficient depth and diameter.

Step 3 - Implant Placement

After creation of the bone bed, install the magnetic implant using the specified screws such that the implant is oriented with the "This Side Up" label facing out. In addition, be sure to place the implant so that the small fiducial mark is located in the most anterior and superior location. Place one screw in each of the bracket holes to securely attach the implant to the skull. The surface of the implant should be nearly flush with the surface of the mastoid. See Figure 3.

Step 4 - Tissue Thinning and Closing the Incision

If the tissue flap is thicker than 6 mm, thin down a 10 mm wide strip of tissue to 5 mm directly over the implant magnets. Suture the incision closed.

Implant Surgical Removal Steps**Step 1 - Making The Incision**

The implant is 6 cm from the ear canal at a 45-degree angle posterior and superior to the ear canal. The incision should be made 7.5 to 8 cm from the external ear canal (1.5 to 2 cm from the location of the implant so the incision is not directly over the implant). After making the incision, expose the skull where the implant has been placed. The incision must be large enough to allow removal of the implant. See Figure 1 above.

Step 2 - Implant Removal

Remove each of the 5 screws holding the magnetic implant in place. It may be necessary to excise bone or tissue growth occurring over the face of the implant. Once the screws are removed the implant can be lifted carefully from the bone or bone beds. Suture incision closed.

Alpha (M) System: Instructions for Placement / Sizing of the Magnetic Spacers and Alpha Sound Processor**Magnetic Spacer Selection**

As with the Headband and Softband, application of appropriate pressure to the mastoid by the Sound Processor is essential for proper functioning of the system and patient comfort. The Alpha (M) configuration is supplied with five strengths of magnetic spacers to allow the fitting physician to adjust the force applied to the mastoid.

The spacers are labeled with "dots" to indicate increasing magnetic strength. Magnetic strength increases based on the number of dots present on the spacer. For example, one dot is a lower strength and five dots is a higher strength.

Several spacer magnet strengths may be applicable for an individual patient. The patient should be prescribed the lowest dot magnet strength that is both comfortable, functions correctly and fits with their age and lifestyle. To determine the appropriate magnetic spacer for a patient use the Magnetic Retention Scale. The Instructions for Use of the scale are included in the scale packaging.

Warning

- Excessive, long term use of magnetic spacers that are too strong may lead to pain and or pressure necrosis for the patient. Always recommend to the patient the lowest magnetic strength spacer that provides sufficient retention.

Caution

- The patient should be instructed to select a lower strength spacer if the pressure on the skin is uncomfortable or if pain, numbness, or tingling occurs. If discomfort continues, advise the patient to discontinue wear and visit the clinic for evaluation and refitting.

Placement of the Alpha (M) Audio Processor

Approximately 4 weeks after surgery, after the incision site has healed and all signs of inflammation are gone, the patient may be fitted with the Alpha (M) Audio Processor.

Warning

- Do not use of the Alpha (M) configuration of the device if the skin contact area of the sound processor is compromised in any way (i.e. – unhealed incision, psoriasis, or eczema). Failure to follow this warning may result in prolonged healing, further irritation, or infection.

4 TECHNICAL DATA

Information about Magnets

Magnets in the Alpha Headband, Softband, and Magnetic Implant create a permanent magnetic field. A permanent magnetic field is not the same as, and should not be compared to, an electromagnetic field such as that of a mobile telephone or a high-voltage transmission line.

Warning

- Both the Alpha (S) and Alpha (M) configurations of the device contain magnets which could impact the performance of magnetically sensitive medical devices such as a pacemaker, an implantable cardioverter defibrillator (ICD), or a ventricular shunt. The magnets of the softband/headband of the Alpha (S) and the magnetic spacer of the Alpha (M) should therefore be kept away from such devices.

Caution

- Keep Otomag™ System components away from magnetic data storage media and electronic devices. The Alpha Sound Processor, Magnetic Implant, Headband, and Softband all contain magnets that can potentially damage magnetic data storage media and electronic devices.
- To avoid potential of damage to the device, do not store the Otomag™ System near magnets or magnetic devices other than the magnetic spacer and implant.
- Patients with a pacemaker or ICD should consult their care provider for safe use of the Otomag system with their cardiac medical device.

Patient-contacting Materials

- Implant Material: The implant consists of a Samarium Cobalt magnet sealed in a biocompatible titanium case.
- Softband Material: The Softband is constructed of an elastic band and consists of approximately 85% nylon/15% Spandex.
- Portion of Headband or Softband that presses against the mastoid (bony protrusion on the skull located directly behind the ear): Fotoplast Lacquer or ABS plastic.
- Implant Template: The implant template is composed of anodized aluminum.
- Surgical Template: The surgical template is composed of polyethylene plastic.
- NO components or packaging for the Otomag™ Bone Conduction Hearing System contain natural rubber latex.

Warning


- The implant is provided clean and non-pyrogenic, but is not sterile. The implant is intended to be sterilized via steam sterilization at the healthcare facility immediately prior to implantation.

Alpha Sound Processor Technical Characteristics

- Magnetically coupled bone conduction sound processor
- Fully programmable digital hearing system
- Programs (programmable by audiologist)
- 8 channels
- 16 frequency bands
- Low battery warning
- Automatic noise reduction
- Automatic feedback cancellation
- Volume control
- Outer material of Alpha (S): plastic with a corrosion resistant steel disk
- Weight: 11 g (without battery) (Alpha 2)


Alpha Sound Processor Performance Data

Measured according to DIN IEC 60118-9 Hearing Aids Part 9: Methods of measurement of characteristics of hearing aids with bone vibrator output.

OFL (Output force level) at 90 dB SPL – Top Line
OFL (Output force level) at 60 dB SPL – Bottom Line

Measured according to DIN IEC 60118-9. Averages include all configurations of softband, headband, and magnetic spacers for a given normal force.

- Acousto – mechanical gain at 60 dB SPL and 1600 Hz: 29 dB
- Maximum gain at 60 dB SPL: 38 dB re 1 µN
- Peak Output at 90 dB SPL: 115 dB re 1 µN
- Frequency range per ANSI S3.22-2009: 280-5000 Hz
- Frequency Processing Range: 125-8000 Hz
- Total Harmonic Distortion per ANSI S3.22-2009: <3%
- Equivalent input noise: 17 dB SPL
- Battery current drain: 0.95 mA
- Processor: 16 band, 8 Channel WDRC
- Crossover Frequency Options: 0-1.25, 0.25-2.75, 1.75-7.75 kHz
- Compression Attack Time: 0.25 ms
- Compression Release Time: 8192 ms
- Patient Volume Control Wheel: 20 dB range
- Power Supply: Alpha 1 – 1.3V, #675; Alpha 2 – 1.3V, #13 (Zinc Air Battery)

5 STORAGE AND HANDLING

Storage

Caution

- To avoid damage to the device, do not subject the Alpha Sound Processor to temperatures below -20° C (-4° F) or above +50° C (122° F) at any time.
- To avoid potential of damage to the device, do not store the Otomag™ System near magnets or magnetic devices other than the magnetic spacer and implant.

Handling

Caution

- Prior to use, inspect packaging to verify that it is intact and has not been damaged. If the packaging is damaged, do not use the device and return it to Sophono. Please see the address and phone number for Sophono at the end of this Physician Manual.

6 PROFESSIONAL SUPPORT

For more information or to report concerns, device malfunctions, adverse events, or patient complaints, please contact:

Sophono, Inc.
5744 Central Ave #100
Boulder, Colorado 80301
USA
Tel: +1 720 407 5160
Fax: +1 720 407 5168
info@sophono.com

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©Copyright Sophono, Inc.
MDSS GmbH